

Changing the battery

When the battery needs changing the hearing aid will start making bleeping sounds. If the hearing aid is not functioning, try changing the battery even if you haven't heard any bleeps.

1. **Open the battery door and remove the old battery.**



2. **Take a new battery from the pack and remove the sticker. Batteries are harmful if ingested. Please carefully store batteries and hearing aids out of reach of small children and animals.**



3. The battery has a flat side with a '+' engraved on top, and a rim on the bottom.



Correct

Place the battery into the battery door from above. The rim side of the battery should be face-down.



Incorrect

Do not force the battery in from the side of the battery door, as this will break the battery door.



Fully closed



Not fully closed



Ensure the battery door is fully closed.

If it is very difficult to close the battery door, recheck the battery is sitting securely.

Accessibility

To view this information in a different language or use text-to-speech reader visit **www.uhb.nhs.uk**, click the yellow and black circular icon in the bottom right of the web page to open the ReachDeck toolbar and then use the search bar to search by the name of the leaflet. If you require this information in another format such as braille, please email **interpreting.service@uhb.nhs.uk**.



How did we do? 😊 😐 😞

If you have recently used our services we'd love to hear about your experience. Please scan the QR code or follow the link to share your feedback to help us improve our services. **Thank you. www.uhb.nhs.uk/fft**

